

## Norco - Oedema Glove - Over Wrist



- Made of nylon / lycra with seams on the outside for smooth oedema control
- Mild compression 15-25 mmHg
- Closed finger style available in wrist length over the wrist style
- Open finger styles available in over-the-wrist style
- Sold individually

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information

Code	Display Name
FBS301151	Norco Oedema Full Finger (Over Wrist) - Left - X Sml
FBS301152	Norco Oedema Full Finger (Over Wrist) - Right - X Sml
FBS301153	Norco Oedema Full Finger (Over Wrist) - Left - Sml
FBS301154	Norco Oedema Full Finger (Over Wrist) - Right - Sml
FBS301155	Norco Oedema Full Finger (Over Wrist) - Left - Med
FBS301156	Norco Oedema Full Finger (Over Wrist) - Right - Med
FBS301157	Norco Oedema Full Finger (Over Wrist) - Left - Lge

Code	Display Name
FBS301158	Norco Oedema Full Finger (Over Wrist) - Right - Lge
FBS301162	Norco Oedema Open Finger (Over Wrist) - Left - X Sml
FBS301163	Norco Oedema Open Finger (Over Wrist) - Right - X Sml
FBS301164	Norco Oedema Open Finger (Over Wrist) - Left - Sml
FBS301165	Norco Oedema Open Finger (Over Wrist) - Right - Sml
FBS301166	Norco Oedema Open Finger (Over Wrist) - Left - Med
FBS301167	Norco Oedema Open Finger (Over Wrist) - Right - Med
FBS301168	Norco Oedema Open Finger (Over Wrist) - Left - Lge
FBS301169	Norco Oedema Open Finger (Over Wrist) - Right - Lge