

## JAY Cushion - LowZone Gel Switch



The any way up gel cushion that offers pressure reduction using gel below a visco elastic layer of foam.

- **Foam & Gel Working In Sync** - The fully enclosed gel pouch benefits from foam and gel working in sync to reduce pressure. Slow moving gel consistency prevents migration and provides comfort and pressure relief. The high performance visco elastic layer of foam is temperature sensitive, encloses the core of the cushion and offers great pressure reduction across the cushion, offering all day support. Orientation marking on cushion as well as cover ensure the correct placement.
- **Wipedown Cover** - An extra protective wipedown cover is easy to care for as it can be wiped clean and washed at 60 degrees while inner Platition cover protects the foam.
- **Pressure Reduction Gel Cushion** - Two-sided cushion can be flipped and rotated to wear evenly. Memory Foam top and bottom layer encasing central gel pouch.

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC536082	JAY Cushion - LowZone Gel Switch - 17 x 17in	17 in	17 in	3 in	158 kg
PTC536084	JAY Cushion - LowZone Gel Switch - 18 x 16in	16 in	18 in	3 in	158 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC536086	JAY Cushion - LowZone Gel Switch - 20 x 18in	18 in	20 in	3 in	158 kg
PTC536088	JAY Cushion - LowZone Gel Switch - 22 x 18in	18 in	22 in	3 in	158 kg