

## Economy Clean Linen Trolley - Closed



- Mobile clean linen trolley
- Transport clean linen in the hospital or aged care environment
- Lightweight construction
- 125 mm castors, wall buffer discs, two shelves, top tray and push handle



### Compliance

No applicable compliance standard

## Warranty Information

- 3 Years

Code	Display Name	Total Length	Total Width	Total Height
LAC353620	Clean Linen Trolley - Small	920 mm	550 mm	1040 mm
LAC353590	Clean Linen Trolley - Large	1100 mm	550 mm	1400 mm
LAC353630	Front Cover - Small	-	-	-
LAC353600	Front Cover - Large	-	-	-