

4H227 Economy Clean Linen Trolley - Closed



- Mobile clean linen trolley
- Transport clean linen in the hospital or aged care environment
- Lightweight construction
- 125 mm castors, wall buffer discs, two shelves, top tray and push handle



Compliance

No applicable compliance standard

Warranty Information

- 3 Years



Code	Display Name	Total Length	Total Width	Total Height	Material
LAC353620	4H227S Clean Linen Trolley - Small	920 mm	550 mm	1040 mm	Aluminium, Rubber
LAC353590	4H227L Clean Linen Trolley - Large	1100 mm	550 mm	1400 mm	Aluminium, Rubber
LAC353630	Front Cover - Small	-	-	-	-
LAC353644	Front Cover - Small	-	-	-	-
LAC353600	Front Cover - Large	-	-	-	-