

## DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT



- Contoured cotton/elastic construction provides comfortable compression.
- Circumferential strap for improved wrist stabilisation.
- Preformed aluminium stay provides anatomically correct fit.
- Contact closure. Application:
- Support following Mild Sprains and strains & cast removal
- Ideal for carpal tunnel syndrome and immobilization of the wrist

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS305380	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - RIGHT - XS
FBS305390	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - RIGHT - SML
FBS305400	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - RIGHT - MED
FBS305410	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - RIGHT - LGE
FBS305420	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - RIGHT - XL
FBS305430	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - LEFT - XS
FBS305440	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - LEFT - SML

Code	Display Name
FBS305450	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - LEFT - MED
FBS305460	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - LEFT - LGE
FBS305470	DONJOY PROCARE 79-8728/9 COMFORTFORM WRIST SUPPORT - LEFT - XL