

Lightweight Cutlery



- These stainless steel utensils are mounted on plastic stems and sealed into light grey foam handles, to give help for those with a weak or painful grip
- The foam is comfortable to hold and washable
- Handle Length: 115 mm
- Handle Diameter: 28 mm
- Weight: 30 - 45 g

Warranty Information

- 1 Year

- ✔ Contracted DVA product
- ✔ Contract approved

Support @ Home

Function: Eating and drinking

Category: Low risk

Click [here](#) for more information

NDIS

Item Code:

05_150900319_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Unit Weight
DLK276220	Cutlery Set	-
DLK276200	Lightweight Spoon	45 g
DLK276190	Lightweight Fork	34 g

Code	Display Name	Unit Weight
DLK276180	Lightweight Knife	45 g
DLK276210	Lightweight Teaspoon	30 g