

## BodyAssist Arm Sling



- Adjustable buckle closure
- Lightweight fabric
- **Strap & Pouch** - facilitate arm positioning
- **Tube** - Loop design facilitates arm positioning
- Pouch w/ Stabilisation Swathe Straps cross behind shoulders to reduce pressure and increase comfort
- Fits left or right arm

 Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Neck Circumference
FBS284520	BodyAssist Arm Sling - Strap - Small	to 120 mm	-
FBS284530	BodyAssist Arm Sling - Strap - Medium	-	1280 - 1640 mm
FBS284560	BodyAssist Arm Sling - Tube - Medium	-	1280 - 1640 mm
FBS284570	BodyAssist Arm Sling - Tube - Large	-	1640 - 2030 mm
FBS284590	BodyAssist Arm Sling - Pouch - Medium	-	1280 - 1640 mm
FBS284600	BodyAssist Arm Sling - Pouch - Large	-	1640 - 2030 mm
FBS284610	BodyAssist Arm Sling - Pouch - Stabilisation Swathe - Small/Medium	-	-

Code	Display Name	Total Length	Neck Circumference
------	--------------	--------------	--------------------