

Limb Protection - Limbkeepers Full Arm Sleeve



- Non-compression soft & breathable knit
- Latex free
- UPF 50+ Sun Protection
- Anti-microbial - silver infused yarns help curb odour causing bacteria (sleeves only)
- Wicking - channels moisture away from skin
- Seamless - no fasteners
- Stays in place - no adhesives
- Form fitting with 4-way rebound stretch
- Unique blend of performance yarns - keep cool when it's warm/warm when it's cool
- Non-restricting
- Interchangeable design for a perfect fit
- Colourfast - easy care/quick dry
- Soft and breathable knit
- One size fits most
- Unisex
- Sold in pairs Indications:
- A preventative solution to help protect fragile skin on hands from abrasions, skin tears, and bruising
- These non-compression knit gloves provide cushioned comfort without bulk
- Limbkeepers will preserve your skin integrity while giving confidence and dignity to the wearer
- These products help reduce injuries, allow you to continue daily activities, maintain quality of life, and independence

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name	Colour
PTA537450	Limbkeepers - Full Arm Sleeve - Heavy(Regular)	Beige
PTA537460	Limbkeepers - Full Arm Sleeve - Heavy(Regular) - Black	Black
PTA537470	Limbkeepers - Full Arm Sleeve - Heavy(Regular) - White	White
PTA537480	Limbkeepers - Forearm Sleeve - Heavy(Regular)	Beige
PTA537490	Limbkeepers - Forearm Sleeve - Heavy(Regular) - Black	Black
PTA537500	Limbkeepers - Forearm Sleeve - Heavy(Regular) - White	White
PTA537510	Limbkeepers - Forearm Sleeve - Light	Beige
PTA537520	Limbkeepers - Forearm Sleeve - Light - Black	Black
PTA537530	Limbkeepers - Forearm Sleeve - Light - White	White
PTA537640	Limbkeepers - Full Arm Sleeve - Light	Beige
PTA537650	Limbkeepers - Full Arm Sleeve - Light - Black	Black
PTA537660	Limbkeepers - Full Arm Sleeve - Light - White	White