

Repose Contur Acute Recline Chair Overlay



- Repose Contur is a pressure redistribution cushion overlay for riser recliner chairs utilising the Repose tried and trusted inflatable air cell and ‘smart valve’ technologies
- Offering pressure area care to the seat, back and leg sections, it secures to the chair easily and safely with straps and thus it can either be fitted to new chairs or be a straightforward retro-fit for users who already have a chair and require a specialist support surface
- Specially formulated extra resistant to chemical disinfectants, especially chlorine
- Adjustable straps for a secure fit on a range of chair sizes
- Repose tried and trusted single air cell and smart valve technologies.
- The zip has an integral cover flap for added protection against fluid ingress
- Provides pressure redistribution across the whole surface for the at-risk areas of the body – sacrum, buttocks, scapulae, vertebrae, shoulders, back of head
- Sophisticated technology in a simple to use recliner chair overlay cushion
- Secure strap system allows for straightforward retro-fit or new chair use
- The cover can be wiped down or machine washed at 60°C
- Cost effective method to provide pressure area care to a recliner chair
- Dimensions provided are when the overlay is inflated

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest	Backrest Height
PTM636654	Repose Contur Acute Recline Chair Overlay	1715 mm	490 mm	50 mm	139 kg	425 mm	490 mm	480 mm	480 mm	745 mm

