

Insulin Cooling Bag - Zip



Glucology cooling pouches are designed to keep your insulin cool and protected for up to 36 hours. Unlike other cooler bags that require batteries or ice packs, this Glucology Pouch only needs to be immersed in cold water to be activated. They can also be reused as many times as you like by re-immersing it in cold water. With its lightweight and compact design, the cooler pouch fits perfectly in your backpack or bag making them perfect when on the move. Available in different sizes and colours.

- Immerse in water to activate and keep your insulin cool for up to 36 hours
- Reusable whenever you need to extend the insulating effects by immersing again in cold water
- Lightweight and compact design makes it travel-friendly and discreet while easily fitting up to 2 syringes and 2 vials at a time
- 12-month warranty, if within 12 months of purchase you are not satisfied with the quality or performance, let us know so we can provide a full refund or product replacement.

✓ Contracted DVA product

