

## Original Highback Sling - Plastic Coated net



- Liko Original HighBack Sling is a basic model which is designed to adapt to the patient without individual adjustments and is one of our most used slings
- It provides for a slightly semi-reclined sitting posture and excellent support for the entire body, which is good for patients with reduced torso stability
- A properly fitted and carefully applied sling will give the patient an excellent sense of security, as well as a high degree of comfort
- The Original HighBack Sling is also recommended for lifting to or from the floor, since it provides a comfortable head support both in sitting and in lying position

 Contract approved

**ndis**

**Item Code:**

05\_123603164\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year



**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394676	Original Highback Sling - Plastic Coated net - Medium Slim	200 kg	Plastic Coated Net