

Vicair Liberty Profile Air Exchange



- Vicair Liberty Profile is an easy to use and comfortable cushion for occasional use
- Designed for comfort in combination with a low to medium level of support and pressure ulcer prevention
- The Vicair Liberty Profile has an integrated foam base with a semi-circle cut out at the location of the sitting bones
- Vicair Liberty cushion with more stability thanks to the shape of the integrated foam base
- Air valves automatically adapt the cushion to your body shape and weight
- Two compartment layout limits sliding
- Easy to clean and disinfect
- Comfortable and rain proof
- An occasional use wheelchair cushion
- For extra comfort the Vicair Liberty Profile is available with an Air Exchange cover. This cover is made of 3D mesh fabric, to make it extra breathable and comfortable
- Incotec Cover available for use in all situations where incontinence and hygiene is an issue, fluid and incontinence resistant, machine washable



Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC609360	Vicair Liberty Profile Cushion - Air Exchange - 18x16in	16 in	18 in	6 cm	150 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
PTC609370	Vicair Liberty Profile Cushion - Air Exchange - 18x18in	18 in	18 in	6 cm	150 kg
PTC609380	Vicair Liberty Profile Cushion - Air Exchange - 18x20in	20 in	18 in	6 cm	150 kg
PTC609390	Vicair Liberty Profile Cushion - Air Exchange - 20x18in	18 in	20 in	6 cm	150 kg
PTC609400	Vicair Liberty Profile Cushion - Air Exchange - 20x20in	20 in	20 in	6 cm	150 kg
PTC609410	Vicair Liberty Profile Cushion - Incotec Cover - 18x16in	16 in	18 in	6 cm	150 kg
PTC609420	Vicair Liberty Profile Cushion - Incotec Cover - 18x18in	18 in	18 in	6 cm	150 kg
PTC609430	Vicair Liberty Profile Cushion - Incotec Cover - 18x20in	20 in	18 in	6 cm	150 kg
PTC609440	Vicair Liberty Profile Cushion - Incotec Cover - 20x18in	18 in	20 in	6 cm	150 kg
PTC609450	Vicair Liberty Profile Cushion - Incotec Cover - 20x20in	20 in	20 in	6 cm	150 kg