

Liko SafetyVest (Mod. 93)



- Sit-to-stand sling designed to assist patients with poor stability
- Self-regulating tensioning function ensures both caregiver and patient feel secure
- Padded fabric provides a high degree of friction and effectively prevents gliding
- Optional Sheepskin padding is available for even more comfort
- Available in Polyester

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394160	SafetyVest Mod 93 Sling - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394170	SafetyVest Mod 93 Sling - Polyester - LARGE	200 kg	Polyester (Loops)