

## Meiko Basket



- Basket for small items with hinged lid is designed to be insert in SUH
- Basket insert with hinged lid is designed to be placed in basket for small items

 Contract approved

### Support @ Home

**Function:** Collecting urine and faeces

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Length	Total Width	Total Height
UPA400100	Meiko Basket Insert	280 mm	100 mm	50 mm
UPA400110	Meiko Basket - for small items	310 mm	240 mm	72 mm