

DONJOY PROCARE 81-8110 SILICONE HEEL CUP



- Molded medical-grade soft silicone.
- Softer center in heel region to help reduce strain to area sensitive to pain.
- Sold in pairs. Application:
- Help reduce shock and vibration loads to joints
- Plantar Fasciitis
- Heel Spurs
- Achilles Tendonitis

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS305650	DONJOY PROCARE 81-8110 SILICONE HEEL CUP - SML/MED
FBS305660	DONJOY PROCARE 81-8110 SILICONE HEEL CUP - LGE/XL