

Freedom Bath Transfer Bench



- Removable backrest can be attached to either side of the bench
- Detachable armrest
- Height adjustable legs. Adjustment pins are spaced at 25 mm increments on three extension legs and 12.5 mm increments on the fourth to maximise stability
- Thicker, more durable powder coated aluminium frame to prevent scuffs and scratches
- Heavy duty sure grip suction rubber tips
- Purchased separately: Extension Legs (GSP650100)

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 1 Year

Compliance

AS/NZS 3973:2009



Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTS153100	Freedom Bath Transfer Bench - HBA424	520 mm	785 mm	825 - 975 mm	200 kg	6 kg	410 mm	680 mm	490 - 640 mm	-
BTS153120	Freedom Bath Transfer Bench w/Full Length Seat - HBA426	520 mm	785 mm	825 - 975 mm	160 kg	6.2 kg	410 mm	680 mm	490 - 640 mm	490 - 640 mm