

Freedom Toilet Seat Raiser - HBA444



freedom
HEALTHCARE

- Padded armrests can swing up and out of the way or be easily removed for easy cleaning
- Entire seat flips up so the bowl can be accessed without the seat (except on lowest height setting)
- Height can be adjusted to 5, 10 or 15 cm to ensure hips remain above knee height
- Attaches to existing toilet seat mounting points

✓ Contracted DVA product

✓ Contract approved

Compliance

AS/NZS 3973:2009

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth
BTT164100	Freedom Toilet Seat Raiser	242	1	380 mm	50, 100 or 150 mm	136 kg	3.8 kg	430 mm