

Spex Mantaray Back Support



- A lightweight solid aluminium back support to support minimal deformities or stability needs, with quick-release hardware
- Minimal to moderate postural support, for mild-moderate torso asymmetries
- Greater functionality as there is minimal hardware between the client and the lateral supports
- Slick and non-obtrusive, provides right amount of spine and lateral stability for different postural needs
- Fully adjustable and repositionable lateral 'wings' means that the therapist/clinician can ensure the right support in the right place for the user
- Improved trunk/pelvic positioning whilst maintaining propelling efficiency
- Minimal contouring, and includes positioning kit for further postural enhancements

 Contract approved

Support @ Home

Function: Accessories for wheelchair seating

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 2 Years

Code	Display Name	Safe Working Load	Backrest Height	Backrest Width
PTB535200	Mantaray Back Support For 12in X 16in Complete With Repositionable Laterals	90 kg	16 in	12 in
PTB535202	Mantaray Back Support For 13in X 16in Complete With Repositionable Laterals	90 kg	16 in	13 in
PTB535204	Mantaray Back Support For 14in X 18in Complete With Repositionable Laterals	90 kg	18 in	14 in
PTB535206	Mantaray Back Support For 15in X 18in Complete With Repositionable Laterals	90 kg	18 in	15 in

Code	Display Name	Safe Working Load	Backrest Height	Backrest Width
PTB535208	Mantaray Back Support For 16in X 20in Complete With Repositionable Laterals	160 kg	20 in	16 in
PTB535210	Mantaray Back Support For 17in X 20in Complete With Repositionable Laterals	160 kg	20 in	17 in
PTB535212	Mantaray Back Support For 18in X 20in Complete With Repositionable Laterals	160 kg	20 in	18 in
PTB535214	Mantaray Back Support For 19in X 20in Complete With Repositionable Laterals	160 kg	20 in	19 in
PTB535216	Mantaray Back Support For 20in X 20in Complete With Repositionable Laterals	160 kg	20 in	20 in