

iCare Bed Rail Covers



- iCare padded rail covers for protection on the Icare side rails
- These padded covers are soft and covered with wipeable medical fabric.
- Valcro holds padded section in place

 Contract approved

Support @ Home

Function: Beds and bed equipment

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Total Width	Total Height
BEA016460	iCare Bed Rail Covers - Bed Stick	200 mm	450 mm
BEA016465	iCare Bed Rail Covers - High Side	400 mm	580 mm
BEA016470	iCare Bed Rail Covers - Low Side	730 mm	435 mm
BEA016475	iCare Bed Rail Covers - Full Length	1400 mm	385 mm