

Elevating Leg Rest M2S MWC



ottobock.

- A great range of accessories to fit the Ottobock M2S MWC range to improve functionality and user experience

Warranty Information

- 2 Years



Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Code	Display Name
WCA859401	Elevating Legrest - Right - 17 to 18 in
WCA859402	Elevating Legrest - Left - 17 to 18 in
WCA859403	Elevating Legrest - Right -20 in
WCA859404	Elevating Legrest - Left - 20 in