

Solaris ExoSoft - Lymphoedema - Wrist to Axilla



- This circular-knit arm sleeve provides comfortable gradient compression and moderate containment
- Wide bands at the wrist and axilla help the sleeve stay in place and resist rolling. Created to replace daytime compression bandaging, or serve as an alternative or supplement to elastic compression stockings
- Microfiber materials are comfortable and breathable while still providing containment for mild to moderate edema
- Soft and smooth materials are easier for patients to don

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
FBS308970	Solaris ExoSoft - (20-30mmHg) Wrist to Axilla - Tall - Beige - Silicone Top - Medium	Beige