

Dreamline G3 Foam



- The Dreamline G3 cushion is a redefine gel cushion with mild contouring, designed for an active user with moderate-to-high pressure management needs
- The G3 cushion offers the support and benefit of gel without the maintenance and potential risks of a traditional gel cushion
- 4-way stretch water-resistant incontinence cover and an outer breathable moisture wicking cover
- Unique and patented gel cluster packs for maximum skin protection
- Includes a set of 10 BAC pads for pelvic obliquity
- 60 mm high contouring with reinforced inserts to bolster the pommel and lateral regions
- With foam Ischial well
- Max user weight: 140 kg

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Safe Working Load	Pressure Care Rating
PTC621100	Dreamline G3 Foam - 12in x 12in	12 in	12 in	140 kg	High Risk
PTC621110	Dreamline G3 Foam - 12in x 14in	14 in	12 in	140 kg	High Risk
PTC621120	Dreamline G3 Foam - 14in x 14in	14 in	14 in	140 kg	High Risk
PTC621130	Dreamline G3 Foam - 14in x 16in	16 in	14 in	140 kg	High Risk

Code	Display Name	Total Length	Total Width	Safe Working Load	Pressure Care Rating
PTC621140	Dreamline G3 Foam - 14in x 18in	18 in	14 in	140 kg	High Risk
PTC621150	Dreamline G3 Foam - 16in x 16in	16 in	16 in	140 kg	High Risk
PTC621160	Dreamline G3 Foam - 16in x 18in	18 in	16 in	140 kg	High Risk
PTC621170	Dreamline G3 Foam - 16in x 20in	20 in	16 in	140 kg	High Risk
PTC621180	Dreamline G3 Foam - 18in x 18in	18 in	18 in	140 kg	High Risk
PTC621190	Dreamline G3 Foam - 18in x 20in	20 in	18 in	140 kg	High Risk
PTC621200	Dreamline G3 Foam - 18in x 22in	22 in	18 in	140 kg	High Risk
PTC621210	Dreamline G3 Foam - 20in x 20in	20 in	20 in	140 kg	High Risk