

## Dublin Wingback Chair



- An eye-catching wingback armchair with a modern feel
- Fresh, clean lines that complement classic or contemporary settings
- Comfortable place to read and relax when placed in a lounge area or bedroom
- Quality construction ensures this armchair will withstand daily use
- Warm walnut legs give a refined look to the armchair
- Colour/material to be advised at time of order



### Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Total Height
FUL052150	Dublin Wingback Chair	660mm	790mm	1105mm