

Oppo 1085 Elbow Sleeve



- Provides full support of elbow region in conjunction with reducing of the insertion of the extensor carpi radialis tendon
- Provides support and protection during activities which require extensive elbow motion
- Helps reduce pain from tendonitis and bursitis Application:
- Weak elbow
- Elbow strain/sprain
- Olecranon bursitis

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Elbow Circumference
FBS301321	Oppo - Elbow Sleeve - Neoprene - Sml	232 - 250 mm
FBS3013212	Oppo - Elbow Sleeve - Neoprene - Med	257 - 280 mm
FBS3013213	Oppo - Elbow Sleeve - Neoprene - Lge	283 - 300 mm
FBS3013214	Oppo - Elbow Sleeve - Neoprene - XL	308 - 330 mm