

## Aspire EasySit Cushion



The Aspire EasySit Cushion is designed to enhance everyday seating comfort through innovative dual-layer gel technology, a low-profile design, and user-friendly features. Lightweight, easy to clean, and suitable for a wide range of chairs and mobility seating, it provides practical support, comfort, and independence.



- **Smart Gel Technology:** Dual-layer gel distributes weight evenly to enhance long-term comfort for chairs, wheelchairs, recliners, and other items.
- **Designed for Better Posture and Independence:** Low-profile design maintains natural posture and has minimal impact on seat height, supporting easier transfers and compatibility with wheelchairs, recliners, office chairs and more.
- **Breathable and Easy to Clean:** Fitted with a stretch fabric cover that promotes airflow and helps prevent heat and moisture build-up. The removable cover is machine-washable for convenient cleaning and improved hygiene.
- **Secure Anti-Slip Base:** Textured, rubber-like base grips seating surfaces to help prevent unwanted cushion movement during transfers or repositioning, enhancing safety and stability.
- **Lightweight and Travel-Friendly:** Compact, portable design makes it easy to move between your wheelchair, car, office chair or recliner — perfect for active and mobile users.
- **Durable Construction:** Cushion made from TPE with a breathable polyester cover, designed for long-lasting performance in everyday use.
- **Simple Care Instructions:** Cover is machine-washable; cushion can be cleaned with mild soap solution. Bleach, harsh detergents, heat drying and autoclave cleaning are not recommended.
- **Available Sizes:** Medium (430 × 415 × 25 mm) and Large (480 × 420 × 25 mm) to suit a variety of seating needs.

## Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height
PTC530050	Aspire EasySit Cushion - Medium	415mm	430mm	25mm
PTC530055	Aspire EasySit Cushion - Large	420mm	480mm	25mm