

## ORTHOBRACE OBSPI1 ORTHOSPICA



- Short design for specific thumb immobilisation
- Removable and adjustable aluminium stays
- Immobilization in three planes
- Soft inner lining Application:
- Stable scaphoid fractures
- CMC injury
- Thumb sprains and strains
- Severe De Quervain's Tenosynovitis

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS308200	ORTHOBRACE OBSPI1 ORTHOSPICA - LEFT - SML
FBS308210	ORTHOBRACE OBSPI1 ORTHOSPICA - RIGHT - SML
FBS308220	ORTHOBRACE OBSPI1 ORTHOSPICA - LEFT - MED
FBS308230	ORTHOBRACE OBSPI1 ORTHOSPICA - RIGHT - MED
FBS308240	ORTHOBRACE OBSPI1 ORTHOSPICA - LEFT - LGE
FBS308250	ORTHOBRACE OBSPI1 ORTHOSPICA - RIGHT - LGE