

## Dreamline STX Gel



- The Dreamline STX cushion is a fully adjustable cushion with full height contouring, designed for inactive users who have asymmetrical positioning requirements and high pressure management needs
- The STX cushion provides the end user the highest amount of adjustability without compromise
- Includes a set of 10 BAC pads for pelvic obliquity
- 4-way stretch water-resistant incontinence cover and an outer breathable moisture wicking cover
- With gel Ischial well
- 80 mm high contouring
- Adjustable gel/foam cluster packs
- Safe working load: 140 kg

 Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Safe Working Load	Pressure Care Rating
PTC621700	Dreamline STX Gel - 12in x 12in	12 in	12 in	140 kg	Very High Risk
PTC621710	Dreamline STX Gel - 12in x 14in	14 in	12 in	140 kg	Very High Risk
PTC621720	Dreamline STX Gel - 14in x 14in	14 in	14 in	140 kg	Very High Risk
PTC621730	Dreamline STX Gel - 14in x 16in	16 in	14 in	140 kg	Very High Risk

Code	Display Name	Total Length	Total Width	Safe Working Load	Pressure Care Rating
PTC621740	Dreamline STX Gel - 14in x 18in	18 in	14 in	140 kg	Very High Risk
PTC621750	Dreamline STX Gel - 16in x 16in	16 in	16 in	140 kg	Very High Risk
PTC621760	Dreamline STX Gel - 16in x 18in	18 in	16 in	140 kg	Very High Risk
PTC621770	Dreamline STX Gel - 16in x 20in	20 in	16 in	140 kg	Very High Risk
PTC621780	Dreamline STX Gel - 18in x 18in	18 in	18 in	140 kg	Very High Risk
PTC621790	Dreamline STX Gel - 18in x 20in	20 in	18 in	140 kg	Very High Risk
PTC621800	Dreamline STX Gel - 18in x 22in	22 in	18 in	140 kg	Very High Risk
PTC621810	Dreamline STX Gel - 20in x 20in	20 in	20 in	140 kg	Very High Risk