

IBICI 'Segreta 140' Womens Compression Stockings



IBICI Segreta 140 (for Women)

- Fashionably sheer, in a range of colours and styles, including pantyhose designed for the shorter, fuller figure
- Maternity pantyhose are not on the DVA/RAP schedule

IBICI general overview:

Designed and manufactured in Italy, IBICI's compression legwear range is ideal for mild vascular conditions, or as a preventative measure. All IBICI garments have a closed toe (full foot). The range consists of:

- Segreta 140 (for Women)
- Segreta Plus (for Women)
- Repomen (for Men)
- Skywalk (Unisex)

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS641879	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 1) - Bge	Beige
PTS641879BE2	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 2) - Bge	Beige
PTS641879BE3	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 3) - Bge	Beige
PTS641879BE4	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 4) - Bge	Beige
PTS641879BE5	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 5) - Bge	Beige
PTS641879BE6	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 6) - Bge	Beige

Code	Display Name	Colour
PTS641879BL1	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 1) - Blk	Black
PTS641879BL2	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 2) - Blk	Black
PTS641879BL3	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 3) - Blk	Black
PTS641879BL4	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 4) - Blk	Black
PTS641879BL5	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 5) - Blk	Black
PTS641879BL6	Ibici Segreta 140 - Knee - 15-20mmHg - C/Toe - (Sz 6) - Blk	Black
PTS641882	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 1) - Bge	Beige
PTS641882BE2	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 2) - Bge	Beige
PTS641882BE3	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 3) - Bge	Beige
PTS641882BE4	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 4) - Bge	Beige
PTS641882BE5	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 5) - Bge	Beige
PTS641882BL1	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 1) - Blk	Black
PTS641882BL2	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 2) - Blk	Black
PTS641882BL3	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 3) - Blk	Black
PTS641882BL4	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 4) - Blk	Black
PTS641882BL5	Ibici Segreta 140 - Waist (regular) - 15-20mmHg - P/hose - C/Toe - (Sz 5) - Blk	Black
PTS641883BE2	Ibici Segreta 140 - Waist (wider) - 15-20mmHg - P/hose - C/Toe - (Sz 2) - Bge	Beige

Code	Display Name	Colour
PTS641883BE3	Ibici Segreta 140 - Waist (wider) - 15-20mmHg - P/hose - C/Toe - (Sz 3) - Bge	Beige
PTS641883BE4	Ibici Segreta 140 - Waist (wider) - 15-20mmHg - P/hose - C/Toe - (Sz 4) - Bge	Beige