

Electric 4HC Bariatric Tilt in Space Shower Commode



- Robust design, manufactured from stainless steel
- Multi height vertical push handle
- Handmade clip-on pressure reducing backrest, removable for cleaning and improved infection control
- Polyurethane swings up armrests
- Generous 700 mm seat width
- High 250kg SWL, ideal for the bariatric environment
- Adjustability through the backrest and seat to suit changing user needs
- Multi position headrest for improved comfort during TIS positioning
- Removable rechargeable battery pack
- Water resistant IPX6 electric actuators for a smooth and effortless tilt action
- Handset controller for freedom to operate around commode

 Contract approved

Support @ Home

Function: Washing, bathing and showering

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 5 Years

Code	Display Name	Total Width	Total Height	Safe Working Load	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BTC069507	Electric 4HC Bariatric Tilt in Space Shower Commode - Attendant Propelled - Swing Away Footrest	700 mm	1415 mm	250 kg	480 mm	700 mm	590 mm	590 mm