

DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT



- Contoured cotton/elastic construction provides comfortable compression.
- Circumferential strap for improved wrist stabilisation.
- Preformed aluminium stay provides anatomically correct fit.
- Contact closure. Application:
- Support following Mild Sprains and strains & cast removal
- Ideal for carpal tunnel syndrome and immobilization of the wrist

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Wrist Circumference
FBS305280	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - XS	11.5 - 14 cm
FBS305290	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - SML	14 - 16.5 cm
FBS305300	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - MED	16.5 - 19 cm
FBS305310	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - LGE	19 - 21.5 cm
FBS305320	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - XL	21.5+ cm
FBS305330	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - XS	11.5 - 14 cm
FBS305340	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - SML	14 - 16.5 cm

Code	Display Name	Wrist Circumference
FBS305350	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - MED	16.5 - 19 cm
FBS305360	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - LGE	19 - 21.5 cm
FBS305370	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - XL	21.5+ cm