

## DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT



- Contoured cotton/elastic construction provides comfortable compression.
- Circumferential strap for improved wrist stabilisation.
- Preformed aluminium stay provides anatomically correct fit.
- Contact closure. Application:
- Support following Mild Sprains and strains & cast removal
- Ideal for carpal tunnel syndrome and immobilization of the wrist

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name
FBS305280	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - XS
FBS305290	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - SML
FBS305300	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - MED
FBS305310	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - LGE
FBS305320	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - RIGHT - XL
FBS305330	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - XS
FBS305340	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - SML

Code	Display Name
FBS305350	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - MED
FBS305360	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - LGE
FBS305370	DONJOY PROCARE 79-8715/6 CTS WRIST SUPPORT - LEFT - XL