

## Clean Linen Trolley



- Mobile clean linen trolley
- Transport clean linen in the hospital or aged care environment
- Lightweight aluminium body and shelves
- 200 mm rubber castors, full wrap wall bumper, three shelves, top multi height push handle
- Pictured with drive unit (optional)



### Compliance

No applicable compliance standard

## Warranty Information

- 3 Years

Code	Display Name	Total Length	Total Width	Total Height
LAC353560	Clean Linen Trolley	1290 mm	640 mm	1440 mm
LAC353570	Clean Linen Trolley Front Cover - Large	-	-	-
LAC353575	Clean Linen Trolley Tow Hitch - Large	-	-	-