

## Oppo Hand/Wrist Positioning Splint



- Two unstretchable straps give greater wrist stability
- Incorporating a molded metal palmer stay which is used to position the wrist and hand

 Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

| Code        | Display Name  | Total Length |
|-------------|---|--------------|
| FBS301470   | Oppo Hand/Wrist Positioning Splint - Right - Small  | 130 - 170 mm |
| FBS301470LL | Oppo Hand/Wrist Positioning Splint - Left - Large   | 216 - 305 mm |
| FBS301470LM | Oppo Hand/Wrist Positioning Splint - Left - Medium  | 170 - 216 mm |
| FBS301470LS | Oppo Hand/Wrist Positioning Splint - Left Small     | 130 - 170 mm |
| FBS301470RL | Oppo Hand/Wrist Positioning Splint - Right - Large  | 216 - 305 mm |
| FBS301470RM | Oppo Hand/Wrist Positioning Splint - Right - Medium | 170 - 216 mm |