

## iCare Mattress - IC15



- The IC15 mattress features the ActiveX™ 84 material which is the main comfort layer in this mattress. ActiveX™ is a high grade of elastic type foam which is responsive to temperature around inflamed areas. ActiveX™ will become softer as it becomes warmer. So around a pressure point it will always soften and allow the blood to circulate freely
- XP Support foam is a high-density foam that is used in the bottom layer of the IC15. XP Support foam is an open cell material that breaths to allow air circulation but is high density to prevent bottoming out. This material will easily contour to any profiling bed
- The IC15 cover is made with the new Bluetec cooling fabric. This is all way stretch fabric to conform with body and mattress contour. It is under woven with a cooling yarn to regulate body temperature
- Icare mattresses don't have any springs which makes them suitable for any profiling or static mattress platforms
- Firm feel makes it easier for turning
- Pair with an Icare bed or suitable on any base platform
- No turn design reducing manual handling
- Anticrease XP base to allow easy contour for any profiling bed

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 10 Year Limited Warranty

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM469729	Mattress IC15 - King Single	2030 mm	1070 mm	150 mm	200 kg

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load
BEM470700	IC15 Mattress - Split Queen	2030 mm	740mm (per side) -		300 kg