

Dreamline G3 Gel



- The Dreamline G3 cushion is a redefine gel cushion with mild contouring, designed for an active user with moderate-to-high pressure management needs
- The G3 cushion offers the support and benefit of gel without the maintenance and potential risks of a traditional gel cushion
- 4-way stretch water-resistant incontinence cover and an outer breathable moisture wicking cover
- Unique and patented gel cluster packs for maximum skin protection
- Includes a set of 10 BAC pads for pelvic obliquity
- 60 mm high contouring with reinforced inserts to bolster the pommel and lateral regions
- With gel Ischial well
- Max user weight: 140 kg

✔ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Safe Working Load	Pressure Care Rating
PTC621300	Dreamline G3 Gel - 12in x 12in	12 in	12 in	140 kg	High Risk
PTC621310	Dreamline G3 Gel - 12in x 14in	14 in	12 in	140 kg	High Risk
PTC621320	Dreamline G3 Gel - 14in x 14in	14 in	14 in	140 kg	High Risk
PTC621330	Dreamline G3 Gel - 14in x 16in	16 in	14 in	140 kg	High Risk

Code	Display Name	Total Length	Total Width	Safe Working Load	Pressure Care Rating
PTC621340	Dreamline G3 Gel - 14in x 18in	18 in	14 in	140 kg	High Risk
PTC621350	Dreamline G3 Gel - 16in x 16in	16 in	16 in	140 kg	High Risk
PTC621360	Dreamline G3 Gel - 16in x 18in	18 in	16 in	140 kg	High Risk
PTC621370	Dreamline G3 Gel - 16in x 20in	20 in	16 in	140 kg	High Risk
PTC621380	Dreamline G3 Gel - 18in x 18in	18 in	18 in	140 kg	High Risk
PTC621390	Dreamline G3 Gel - 18in x 20in	20 in	18 in	140 kg	High Risk
PTC621400	Dreamline G3 Gel - 18in x 22in	22 in	18 in	140 kg	High Risk
PTC621410	Dreamline G3 Gel - 20in x 20in	20 in	20 in	140 kg	High Risk