

Medifab Flux Walker

The Flux is a lightweight but sturdy posterior walking aid which encourages upright posture and promotes natural walking movements.



- Height-adjustable grip bars grow with the user, making it easy to adjust to different body heights
- Folding mechanism with safety lock eases transport and storage
- Friction brakes and bumper wheels included in package
- Opens out in walking direction and optically and physically free from interfering structural parts
- Reverse-roll locks (included in package) can be deactivated so the walker can roll forward and backward
- Grip Height: 510 mm - 880 mm

Code	Display Name	Total Length	Total Width	Safe Working Load	Unit Weight
PED470076	Medifab Flux Walker (as per script)	810 mm	690 mm	75 kg	6.8 kg