

Aspire Mattress Cover



- Modular foam sections can be removed for cleaning and transfer
- Foam side sections conform to profiling beds with head raises / knee breaks
- Heat sealed components for fluid impermeability and maximum infection control
- Vapour permeable, waterproof top cover
- Multi-stretch fabric to minimise friction and reduce shear
- Seams tape sealed for added bacterial protection and infection control

Warranty Information

- 2 Years

 Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

Compliance

IEC 60601-1-2



Code	Display Name	Total Length	Total Width	Total Height
PTM639030	Aspire Concave Positioning Cover 2000x900x125mm	2000 mm	900 mm	125 mm
PTM639035	Aspire Concave Positioning Cover 2000x1060x125mm	2000 mm	1060 mm	125 mm
PTM639040	Aspire Concave Positioning Cover 2000x900x200mm	2000 mm	900 mm	200 mm
PTM639045	Aspire Concave Positioning Cover 2000x1060x200mm	2000 mm	1060 mm	200 mm