

BodyAssist Abdominal Support



- Helps support lowered abdominal lift capacity
- Provides improved body shape and greater sense of confidence when worn
- Provides complete freedom of movement and comfort
- Breathable fabric for maximum all-day comfort and compression
- Simple adjustment with an optimum fit



Code	Display Name	Safe Working Load	Material	Waist/Hip Circumference
FBS283690	BodyAssist Abdominal Support Belt - Medium	220 mm	Elastic	850 - 950 mm
FBS283700	BodyAssist Abdominal Support Belt - Large	220 mm	Elastic	950 - 1050 mm
FBS283710	BodyAssist Abdominal Support Belt - XLarge	260 mm	Elastic	1050 - 1150 mm
FBS280550	BodyAssist Abdominal Support Belt - XXL	260 mm	Elastic	1150 - 270 mm