

Underarm Crutches - Red Dot



- Ergonomic underarm crutches with a patented clip-lock adjustable handle
- Push-button foot pieces, enabling quick, accurate measuring and fitting.
- Lightweight aluminum frame and high-density rubber components ensuring long-lasting wear and exceptional comfort, enhancing your mobility experience.
- Patented clip - lock adjustable handle and push button foot pieces for fast, precise measuring and fitting
- Lightweight aluminium frame
- High density rubber cushions and handgrips providing long wear and greater comfort
- Ferrule: 22 mm
- Product Weight: 1.25 kg
- User Height: WAC692300 (Child) = 1220-1370mm, WAC692500 (Youth) = 1320-1570mm, WAC692200 (Adult Medium) = 1550-1750mm, WAC692400 (Adult Tall) = 1780-1980mm

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_120306087_0105_1_2

Walking Aids -

Composite:05_416_0105_1_2

Click [here](#) for more information

Compliance

ISO 11199-2:2005

Warranty Information

- 1 Year

Code	Display Name	RRP	Box Quantity	Total Height	Safe Working Load
WAC692300	Red Dot Underarm Crutches - Child/Paediatric	53	8	787 - 940 mm	136 kg
WAC692500	Red Dot Underarm Crutches - Youth/Small	77	8	940 - 1140 mm	136 kg
WAC692200	Red Dot Underarm Crutches - Adult/Medium	77	8	1120 - 1320 mm	136 kg

Code		Display Name		RRP	Box Quantity	Total Height	Safe Working Load
WAC692400		Red Dot Underarm Crutches - Tall/Large		92	8	1320 - 1520 mm	136 kg