

Leckey Freestander



- Upright stander designed to assist in standing and weight bearing development
- Highly adjustable to suit individual users
- Includes chest and knee supports, hip belt, angle adjustable tray and activity bowl
- Mobility base with 4 locking castors
- User Standing Height - PED479170 = 800 - 1151 mm, PED479180 = 1100 - 1448 mm, PED479190 = 1400 - 1849 mm
- Chest Pad Height - PED479170 = 470 - 680 mm, PED479180 = 680 - 975 mm, PED479190 = 975 - 1350 mm



ndis

Item Code:

05_053603010_0103_1_2

Click [here](#) for more information

Warranty Information

- 2 Years

LECKEY®

Code	Display Name	Safe Working Load	Unit Weight
PED479170	Leckey Freestander - Size 1, Ages 1.5-5 years	60 kg	26.5 kg
PED479180	Leckey Freestander - Size 2, Ages 4-10 years	80 kg	30.5 kg
PED479190	Leckey Freestander - Size 3, Ages 9-18 years	90 kg	34 kg