

## JAY - Lite



- An ideal off-the-shelf complement to ultra lightweight wheelchairs
- Constructed with lightweight 3DX Microclimatic spacer fabric and Airflow base to provide excellent microclimatic heat and moisture dissipation
- The OptiWell pelvic loading area helps to achieve high levels of pressure distribution and stability in an extremely lightweight cushion
- Available with positioning inserts for those with postural needs (please ask our Aidacare C.A.S.E consultants for help selecting these)
- Available in sizes 14 x 14" to 24 x 24"
- All models include: Incontinence resistance cover

### Warranty Information

- 1 Year Cushion
- 6 Months Cover

# JAY<sup>®</sup>

✓ Contracted DVA product

✓ Contract approved

**ndis**

**Item Code:**

05\_181006190\_0105\_1\_2

**Pressure Care Cushion –**

**Composite:**05\_429\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	Total Height	Safe Working Load	Pressure Care Rating	Material
PTC571913	JAY Lite - 16 x 18in	2.25 in	113 kg	Medium Risk	Spacer Fabric
PTC585157	JAY Lite - 18 x 16in	2.25 in	113 kg	Medium Risk	Spacer Fabric