


Liko Solo Highback Sling (Mod. 25)



- Designed for single patient use
- Can be used in most lifting situations
- Designed to adapt to the patient without any individual adjustments
- Made of strong non-woven polypropylene fabric that is non-washable and engineered to be disposed of when it has become soiled or when the patient no longer requires it
- Supplied in boxes of 20 (Med-Lrg) & 8 (XL)

Warranty Information

- 1 Year



ndis

Item Code:

05_123603164_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394740	Solo Highback Sling - Non-Woven Polypropylene - MEDIUM Slim	200 kg	Polyester (Loops)
LSS394750	Solo Highback Sling - Non-Woven Polypropylene - MEDIUM	200 kg	Polyester (Loops)
LSS394760	Solo Highback Sling - Non-Woven Polypropylene - LARGE	200 kg	Polyester (Loops)
LSS394770	Solo Highback Sling - Non-Woven Polypropylene - XLARGE	500 kg	Polyester (Loops)