

## Liko Solo Highback Sling (Mod. 25)



- Designed for single patient use
- Can be used in most lifting situations
- Designed to adapt to the patient without any individual adjustments
- Made of strong non-woven polypropylene fabric that is non-washable and engineered to be disposed of when it has become soiled or when the patient no longer requires it
- Supplied in boxes of 20 (Med-Lrg) & 8 (XL)

 Contract approved

**ndis**

**Item Code:**

05\_123603164\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year

**Compliance**

AS/NZS ISO 10535:2011



Code	Display Name	Safe Working Load	Material
LSS394740	Solo Highback Sling - Non-Woven Polypropylene - MEDIUM Slim	200 kg	Polyester (Loops)
LSS394750	Solo Highback Sling - Non-Woven Polypropylene - MEDIUM	200 kg	Polyester (Loops)
LSS394760	Solo Highback Sling - Non-Woven Polypropylene - LARGE	200 kg	Polyester (Loops)
LSS394770	Solo Highback Sling - Non-Woven Polypropylene - XLARGE	500 kg	Polyester (Loops)