

ISA XPLUS Stander



- Three foot plate heights available with removable heel block and round shape for toilet approach. Intermediate foot plate supplied with lifter
- The push bar and foot pedal have an ergonomic shape for easy lifter operation
- Sling hooks are flexible and secure for a quick, straightforward sling attachment
- Leg band padded for extra comfort. Easy to exchange and washable
- Curved legs have an ergonomic shape for easy manoeuvring around chairs and toilet bases
- Lockable rear casters are added security and safety for the patient and user
- Removable battery has a second battery charging to reduce lifter down time

 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

| Code | Display Name | Total Length | Total Width | Safe Working Load | Lift Range |
|-----------|--------------------------------------|--------------|---------------|-------------------|--------------|
| LSS378050 | Invacare Lifters - ISA XPLUS Stander | 1280 mm | 640 - 1020 mm | 200 kg | 505 - 675 mm |