

## Aspire MAXI Adjustable Bedside Commode



Experience enhanced convenience and comfort with the Aspire MAXI bedside commode. Crafted for bariatric users with a generously proportioned seat and extensive width adjustment between the arms, it provides ample space for larger users up to 310kg SWL.

- **Superior Durability:** Constructed with high-quality materials, it can withstand substantial weight, making it a reliable choice for MAXI users. The robust frame ensures stability, providing users with a sense of security and peace of mind.
- **Ease of Use:** Soft padded seat is easily lifted for toileting. Flip-up toilet seat for bowl access and cleaning.
- **Adjustable and Versatile:** Height adjustable legs and width adjustable frame means it is easily adjusted to suit a wide range of user's sizes and needs.
- **Bedroom Convenience:** Toileting and a comfortable chair are both conveniently and safely accessible within the bedroom.

 Contract approved

### Support @ Home

**Function:** Toileting

**Category:** Prescribed

Click [here](#) for more information

## Warranty Information

- 1 Year

Code	Display Name	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Colour	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
BEC045125	Aspire MAXI Adjustable Bedside Commode	715 - 915 mm	890 - 1015 mm	310 kg	Vinyl	12.3 kg	Champagne	475 mm	530 mm	490 - 615 mm	490 - 615 mm