

Medi Stocking Applicator



- Simplifies donning of compression stockings
- Especially designed for patients with large stocking sizes
- Can be used with open or closed toe compression stockings
- Designed for stocking wearers with considerably restricted movement
- Extra-wide opening
- Ideal for those with arthritis or limited bending motion.

- ✓ Contracted DVA product
- ✓ Contract approved

Code	Display Name
PTS644920	Stocking Applicator - Medi XXL Wide Calf - wire frame - Size Universal
PTS644930	Stocking Applicator - Medi Variable Handles - wire frame - Size Universal