

## Pressure Relief Chair



- Tilt / Recline / Elevating Leg-rest
- Visco elastic memory foam throughout minimises pressure without compromising postural support
- Careflex three-way stretch fabric minimises sheer (friction) whilst promoting air-flow
- The Anti-Microbial, Anti-Bacterial properties of the fabric help with infection control and skin health
- Tilt-in Space adjustable seat positioning with infinite locking
- Adjustable wings and arms for easy transfer
- Polyurethane, Anti-Microbial, Anti- Bacterial, Three-way stretch fabric
- Adjustable backrest for full recline
- Adjustable legrest with infinite raise and lower positions helps position the user to help minimise fluid retention
- Four low-friction locking castors -Three braking, One directional
- Soft-touch brake and recline levers for easier operation



\*\*\* Make and model of hire items may vary depending on availability. Minimum rental is two weeks.

Code	Display Name
H_CHP198390	Pressure Relief Chair