

ORTHOBRACE OBWRU1-U ORTHOWRIST UNI



- Extra-long forearm design
- 3 straps provide optimal support*
- Volar, dorsal and lateral aluminium stays
- Soft inner lining Application:
- Buckle fractures
- Mild, stable distal radius and wrist fractures
- Support following cast removal
- Wrist sprains and strains

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS308350	ORTHOBRACE OBWRU1-U ORTHOWRIST - UNIVERSAL