

Body Assist PSB Patella Stabilizing Brace



- Designed according to knee cap taping technique, locking the knee cap
- Anatomically shaped to provide orthopaedic support
- Retains some body warmth & can increase blood circulation
- A great brace for anyone who needs to tape their knee cap

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS285280	BodyAssist - Knee Brace - Patella Stabilising Brace - XL	420 mm+
FBS285270	BodyAssist - Knee Brace - Patella Stabilising Brace - Regular	320 - 420 mm