

Cura Alarms - Wrist Transmitter



- 433 MHz frequency
- One-year life Lithium battery (not replaceable)
- On/off using caregiver reset key
- Disposable hospital grade PVC wrist straps
- 4 grams in weight
- Irreversible clasp
- 105 mm length shortest length
- Hospital grade PVC

Warranty Information

- 1 Year



 Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_222704219_0123_1_2

Personal Alarm – Composite:

05_427_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width	Unit Weight
BEA007461	Wristband Transmitter	36 mm	11 mm	0.4 kg